

A cluster of overlapping, semi-transparent geometric shapes in shades of light blue, teal, green, and purple, located in the top-left corner of the slide.

# **ACEMFS FUT Minna Bioinformatics Workshop**

## **Primer Design**

**Itunuoluwa Isewon PhD**  
Covenant University

A cluster of overlapping, semi-transparent geometric shapes in shades of light gray, located in the bottom-right corner of the slide.

A cluster of overlapping, semi-transparent geometric shapes in the top-left corner, including shades of blue, green, red, and purple.

# **Primer Design**

A cluster of overlapping, semi-transparent light gray geometric shapes in the bottom-right corner, resembling a stylized flower or star.

# DNA & PCR Essentials



- **DNA Structure Review:**

- DNA is a double helix composed of two complementary strands.
- Bases: Adenine (A) pairs with Thymine (T); Guanine (G) pairs with Cytosine (C).
- DNA is antiparallel; primers bind 3'→5' on template, extend 5'→3'

- **What is a Primer?**

- A short, synthetic, single-stranded DNA sequence (typically 18-25 bases long).
- They act as starting points for DNA polymerase.
- Designing primers for PCR is an essential wet-lab skill with a bioinformatics foundation.

# DNA & PCR Essentials



- **PCR Basics:**

- **Polymerase Chain Reaction:** An in vitro method to amplify specific DNA segments.
- **Three Steps:** Denaturation, Annealing, and Extension.
- Specificity = primer sequence + annealing temp + template uniqueness
- Amplicon length:
  - Conventional PCR: 100–1000 bp (sweet spot 300–800 bp)
  - qPCR: 70–200 bp (for efficiency)

# Primer Design Rules



- **Length:** 18-25 nucleotides (20–22 nt typical).
- **Melting Temperature ( $T_m$ ):**
  - The temperature at which half of the DNA duplex separates into single strands.
  - **Rule:** Primers should have a  $T_m$  between 55-65°C. For a primer pair, the  $T_m$  should be within a 5°C range of each other.
  - **Calculation:** Varies by formula (e.g., Wallace Rule:  $T_m = 2(A+T) + 4(G+C)$ ). Tools calculate this for you.
- **GC Content:**
  - The percentage of Guanine and Cytosine bases.
  - **Rule:** Aim for 40-60%. High GC content can increase  $T_m$  and lead to stronger binding.

# Primer Design Rules

- **A-T/G-C Clamps:** A G or C at the 3' end strengthens binding. Avoid an A or T.
- **Avoid:**
  - 3 G/Cs at 3' end in a row (non-specific priming)
  - Runs of  $\geq 4$  identical bases
  - Self-complementarity (hairpins) and cross-dimerization ( $\Delta G$  too negative)
- **Amplicon:** pick region without long repeats, extreme GC, or low complexity
- **Check:** secondary structures, off-targets, SNPs, indels

# Melting Temperature (T<sub>m</sub>) Basics



- **Rule-of-thumb (for short primers):**  $T_m \approx 2^{\circ}\text{C} \times (\text{A} + \text{T}) + 4^{\circ}\text{C} \times (\text{G} + \text{C})$
- **Preferred:** nearest-neighbor model with salt/oligo corrections (tool-computed)
- Aim for **annealing**  $\approx T_m - 3$  to  $-5^{\circ}\text{C}$ ; optimize empirically

# Primer Secondary Structures & Specificity



- **Secondary Structures:**

- **Hairpins:** A primer folding back on itself.
- **Self-dimers:** Two identical primers binding to each other.
- **Cross-dimers:** Two different primers binding to each other.
- **Rule:** Avoid all of these! They compete with the target DNA and reduce PCR efficiency.

- **Specificity is Key:**

- Primers must bind only to the intended target sequence.
- **Primer-BLAST:** The essential tool to check for off-target binding in a genome.



# Mycotoxin Pathways & Gene Targets



## Aflatoxins (*Aspergillus flavus/parasiticus*)

- Clustered genes; key regulator **aflR**; PKS **aflC/pksA**; structural genes (**aflD**, **aflM**, ...)
- Targets: **aflR** (regulation), **aflD/aflM** (biosynthesis), **pksA** (polyketide step)

## Trichothecenes (*Fusarium graminearum*, etc.)

- Core **TRI** cluster: **TRI5** (trichodiene synthase), **TRI4**, **TRI6**, **TRI10**
- Targets: **TRI5** for presence/chemotype panels

## Ochratoxin A (*Aspergillus/Penicillium*)

- **ota** cluster (e.g., **otaA/pks**, **otaB**, **otaC**)—species dependent

Use NCBI Nucleotide/Gene to fetch sequences (genomic or mRNA/cDNA)

# Applications & Assay-Specific Considerations

- **End-point PCR (screening/ID):** 300–800 bp; robust primers; gel visualization
- **qPCR/RT-qPCR (quantification):** 80–150 bp; high efficiency (90–110%); MIQE compliance; melt curve or probe
- **Sanger sequencing:** product ~600–900 bp; avoid homopolymers; flank mutation site
- **Mutation/SNP detection:** allele-specific primers or HRM; position variant at/near 3' end for AS-PCR
- **Multiplexing:** balanced  $T_m$ /GC, minimal primer-primer interactions, staggered product sizes

# Primer Design Tools



- **Primer3:**
  - The classic, open-source tool.
  - **Strengths:** Highly customizable, calculates various primer metrics, excellent for designing standard PCR primers.
- **NCBI Primer-BLAST:**
  - Combines primer design with a powerful specificity check against the NCBI database.
  - **Strengths:** Ensures primers are unique and won't amplify off-target sequences.
- **OligoAnalyzer (IDT):**
  - A tool for analyzing existing primers.
  - **Strengths:** Calculates  $T_m$ , checks for hairpins and dimers, provides a comprehensive report on primer quality.
- **Other Tools:**
  - **Geneious:** Integrated software with a user-friendly primer design module.
  - **SnapGene:** Another integrated tool for molecular biology.

# Workflow



- Get target DNA/mRNA sequence (GenBank/RefSeq; confirm species/strain)
- Define region (avoid low complexity; include/exclude sites as needed)
- Design candidate pairs in Primer3 (set product size & T<sub>m</sub>)
- Check specificity in Primer-BLAST (select organism DB; add host crops to “exclude”)
- QC in OligoAnalyzer (hairpins/dimers,  $\Delta G$  thresholds)
- Pick final pair(s); document/export; order
- Bench validation: gradient Ta, efficiency, specificity (melt curve/gel)

# Common Pitfalls & Troubleshooting



## No PCR Product:

- **Potential Cause:** Primers aren't binding (wrong  $T_m$ ), DNA quality is poor, or primers are forming dimers.
- **Troubleshoot:** Re-design primers, check DNA concentration/quality, run a gel to check for dimers.

## Non-Specific Bands:

- **Potential Cause:** Primers are binding to unintended sites in the genome.
- **Troubleshoot:** Increase annealing temperature, re-design primers with higher specificity, use Primer-BLAST.

## Low Yield:

- **Potential Cause:** Sub-optimal primer design, low primer concentration, or poor enzyme performance.
- **Troubleshoot:** Optimize annealing temperature, check primer quality, increase primer concentration.

# Demo 1: Primer3 (Design on Aflatoxin aflR)



**Goal:** Design qPCR-grade primers (100–150 bp) for **aflR** (*Aspergillus flavus*)

**Inputs:** [aflR](#) coding sequence (mRNA or genomic region)

**Constraints:** Length 20–24 nt; Tm  $60 \pm 1$  °C; GC 40–60%; product 100–150 bp

**Deliverables:** 2–3 candidate primer pairs

**Steps:**


**Get sequence:** In NCBI Gene, search “*aflR Aspergillus flavus*”. Open Gene record → pick RefSeq mRNA or genomic sequence; click **FASTA**; copy sequence.

**Open** [Primer3](#)/[Primer3Plus](#) web: paste sequence. Set:

- **Product Size Ranges:** 90-160
- **Primer Size:** 20–24; **Tm:** 59–61 °C; **GC%:** 40–60
- (Optional) Mark **Targets/Included/Excluded Regions** (e.g., avoid low complexity)

Click **Pick Primers**. Review top 5 pairs: lengths, Tm, GC, penalty score.

Copy candidates to a scratch sheet/notepad for downstream checks.



# Primer3web version 4.1.0 - Pick primers from a DNA sequence.

[disclaimer](#)[code](#)[cautions](#)

Select the [Task](#) for primer selection

[Template masking](#) before primer design ([available species](#))

[Select species](#)

[Nucleotides to mask in 5' direction](#)

[Primer failure rate cutoff](#)

[Nucleotides to mask in 3' direction](#)

Paste source sequence below (5'→3', string of ACGTNacgtn -- other letters treated as N -- numbers and blanks ignored). FASTA format ok. Please N-out undesirable sequence (vector, ALUs, LINEs, etc.) or use a [Mispriming Library \(repeat library\)](#)

```
TATCTTGGGATGGGAAGCTGTGGGAGGTAGAGGGTGCCCGAGCTGAAGCATTGGCCGCGACAACTCCG
GGGGCTCGATTCTGGGGCCGCTATCGTGGTGTCCCAACATGACTCTGCTCTCTGGAGGATGG
AAAAATCTCCCCACTGGAATGGATCTTGGAGGGGAAGGAAACGCAGATTGGTTGGTGCGGCACGTCTC
GAATGGAAACCAAACTTGGACTTTGCCGATATGCATAGCCTTGTTCGCCCAAGTCATGGATCTATGAACG
ACGGCCC
```

☒ Pick left primer, or use left primer below

☐ Pick hybridization probe (internal oligo), or use oligo below

☒ Pick right primer, or use right primer below (5' to 3' on opposite strand)

[Pick Primers](#)[Download Settings](#)[Reset Form](#)

[Sequence Id](#)

A string to identify your output.

[Targets](#)

E.g. 50,2 requires primers to surround the 2 bases at positions 50 and 51. Or mark the [source sequence](#) with [ and ]: e.g. ...ATCT[CCCC]TCAT.. means that primers must flank the central CCCC.

[Overlap Junction List](#)

E.g. 27 requires one primer to overlap the junction between positions 27 and 28. Or mark the [source sequence](#) with -: e.g. ...ATCTAC-TGTCAT.. means that primers must overlap the junction between the C and T.

# Primer3Plus

pick primers from a DNA sequence

[More...](#)[Source Code](#)[Help](#)[About](#)

[Load server settings:](#)

**Task:**

Select primer pairs to detect the given template sequence. Optionally targets and included/excluded regions can be specified.

Main

General Settings

Advanced Settings

Internal Oligo

Penalties

Advanced Seq.

[Sequence Id:](#)

[Paste template sequence below](#) or upload sequence file:  No file chosen

```
TATCTTGCGGATGGGAAGCTGTGGGAGGTAGAGGGTGCCCGAGCTGAAGCATTGCCC GCGACAACTCCG
GGGGCTCGATTTCTGGGGCCGCTATCGTGGTGTCCCAACATGACTCTGTCTCTCTCTGGAAGGATGG
AAAATTCTCCCACTGGAAATGGATCTTGGAGGGGAAGGAAACGCAGATTTGGTTGGTGCGGCACGTCTC
GAATGGAACCAAACTTGACTTTGCCGATATGCATAGCCTTGTTGCCCCAAGTCATGGATCTATGAACG
ACGGCCC
```



Mark selected region:

< Excluded >

[ Target ]

{ Included }

Clear

Regions from Seq.

[Excluded Regions:](#) <  >



# Demo 2: Primer-BLAST (Specificity Check)



Goal: Verify **afIR** primer specificity

DB/Organisms:

- **Target include:** *Aspergillus flavus* (taxid 5059)
- **Also test against:** *Aspergillus* spp. (cross-reactivity), host crops (e.g., *Zea mays*, *Triticum aestivum*)

Settings: Default specificity **ON**; choose smallest relevant database (e.g., RefSeq genomic) for precision

Open [NCBI Primer-BLAST](#).

Paste forward/reverse primers (from Primer3) **and** paste the template FASTA (improves alignment context).

Under **Primer Pair Specificity Checking**:

- **Organism:** *Aspergillus flavus*; database: **RefSeq genomic** (or **nr** for broad checks)
- Add additional organisms (maize/wheat/rice) to ensure no host hits

**Get Primers / Submit** → inspect **Amplification products** table, predicted sizes, off-target hits, and alignment details.

**Export:** “Send to → File” (or copy TSV). Keep Job ID for records.



## Primer-BLAST

A tool for finding specific primers

Finding primers specific to your PCR template (using Primer3 and BLAST).

### Primers for target on one template

### Primers common for a group of sequences

#### PCR Template

[Retrieve recent results](#) [Publication](#) [Tips for finding specific primers](#)

[Save search parameters](#)

[Reset page](#)

Enter accession, gi, or FASTA sequence (A refseq record is preferred) ?

[Clear](#)

```
AAAAATTCTCCCACTGGAAATGGATCTTGGAGGGGAAGGAAACGCAGATTGGTTGGTGCGG
CACGTCTC
GAATGGAAACCAAACTTGGACTTTGCCGATATGCATAGCCTTGTTCGCCCAAGTCATGGATCT
ATGAACG
ACGGCCC
```

Or, upload FASTA file

No file chosen

Range ?

[Clear](#)

From

To

Forward primer

Reverse primer

#### Primer Parameters

Use my own forward primer  
(5'->3' on plus strand)

Use my own reverse primer (5'-  
>3' on minus strand)

?

[Clear](#)

?

[Clear](#)

Min

Max

PCR product size

70

1000

# of primers to return

10