



ACEMFS FUT Minna Bioinformatics Workshop

**Sequence Retrieval &
Quality Control using
Galaxy**

Itunuoluwa Isewon PhD
Covenant University



Workshop Outline




By the end of this session, participants will be able to:

- Retrieve sequences of mycotoxins & fungal enzymes from databases.
- Upload & organize datasets in [Galaxy](#).
- Perform QC using Galaxy tools ([FastQC](#))

Introduction to Galaxy



- Web-based platform for bioinformaticians.
 - No coding required.
 - Supports reproducible research.
 - Widely used for genomics, proteomics and transcriptomics analysis.
- 

The image features a white background with decorative geometric elements. In the top-left corner, there is a cluster of overlapping triangles in shades of blue, green, and red. In the bottom-right corner, there is a cluster of overlapping triangles in various shades of gray.

Sequence Retrieval

Sequence Retrieval



NCBI Exploration:

1. Go to [NCBI](#)
2. In the search bar, paste the **lcc9 gene**.
3. Click on the dropdown arrow next to the search bar and select **Nucleotide**.
4. Click the search button.

Sequence Retrieval

NIH National Library of Medicine
National Center for Biotechnology Information

Log in

All Databases ▾ lcc9 gene Search

NCBI Home
Resource List (A-Z)
All Resources
Chemicals & Bioassays
Data & Software
DNA & RNA
Domains & Structures
Genes & Expression
Genetics & Medicine
Genomes & Maps
Homology
Literature
Proteins
Sequence Analysis
Taxonomy
Training & Tutorials
Variation

NCBI
The National Library of Medicine advances science and health by providing access to genomic information.
[Home](#) | [Mission](#) | [Organization](#) | [NCBI News & Blog](#)

Submit
Submit manuscripts and databases

Download
Transfer NCBI data to your computer

Learn
Find help documents, attend a class or watch a tutorial

Develop
Build and code applications

Analyze
Identify an NCBI tool for your data analysis task


Research
Explore NCBI research and collaborative projects

Popular Resources
[PubMed](#)
[Bookshelf](#)
[PubMed Central](#)
[BLAST](#)
[Nucleotide](#)
[Genome](#)
[SNP](#)
[Gene](#)
[Protein](#)
[PubChem](#)


NCBI News & Blog
Upcoming Changes to GenBank Project Lists and Symlinks
20 Aug 2025
The volume of GenBank WGS, TSA, and T1 S data sets is large and continues to grow.
PubMed Central (PMC) Beta Search to Replace Current Search
19 Aug 2025
Coming the week of September 7th The

Sequence Retrieval



5. On the results page click the **first hit**.
 6. What is the name of the **organism**?
 7. How many **base pairs** does it have?
 8. What is its **accession** and **version numbers**?
- 

Sequence Retrieval

 **National Library of Medicine**
National Center for Biotechnology Information

Log in

Nucleotide

Nucleotide

lcc9

Search

Create alert Advanced

Help

Species
Animals (4)
Fungi (21)
Bacteria (1)
Viruses (1)
Customize ...

Molecule types
genomic DNA/RNA (19)
mRNA (8)
Customize ...

Source databases
INSDC (GenBank) (20)
RefSeq (6)
Customize ...

Sequence Type
Nucleotide (27)

Genetic compartments
Mitochondrion (4)
Custom range...

Sequence length
Custom range...

Release date
Custom range...

Revision date
Custom range...

Summary 20 per page Sort by Default order Send to: Filters: [Manage Filters](#)

See [lcc9 laccase-9 precursor](#) in the Gene database
[lcc9](#) reference sequences [Protein](#) (1)

Items: 1 to 20 of 27

<< First < Prev Page 1 of 2 Next > Last >>

☐ [Fusarium oxysporum strain 4287 Lcc9 \(lcc9\) gene, complete cds](#)
1. 3,669 bp linear DNA
Accession: EF990898.1 GI: 152013647
[Protein](#) [Taxonomy](#)
[GenBank](#) [FASTA](#) [Graphics](#)

☐ [Neurospora crassa OR74A Lcc9 \(NCU04593\), partial mRNA](#)
2. 1,635 bp linear mRNA
Accession: XM_952834.3 GI: 758995457
[BioProject](#) [BioSample](#) [Protein](#) [PubMed](#) [Taxonomy](#)
[GenBank](#) [FASTA](#) [Graphics](#)

☐ [Neurospora crassa OR74A linkage group V, whole genome shotgun sequence](#)
3. 6,436,246 bp linear DNA
Accession: NC_026505.1 GI: 759001908
[Assembly](#) [BioProject](#) [BioSample](#) [Protein](#) [PubMed](#) [Taxonomy](#)
[GenBank](#) [FASTA](#) [Graphics](#)

Results by taxon
Top Organisms [\[Tree\]](#)
[Laccaria bicolor](#) (4)
[Neurospora crassa OR74A](#) (3)
[Lentinula edodes](#) (3)
[Agaricus bisporus var. bisporus H97](#) (3)
[Triphophysa robusta](#) (3)
All other taxa (11)
[More...](#)

Find related data
Database: [Select](#)
[Find items](#)

Search details
lcc9[All Fields]
[Search](#)


Activate Windows
Go to Settings to activate Windows.
[See more...](#)

Sequence Retrieval



9. On the right, click on the dropdown arrow next to Send to:
10. Select file under, Choose **Destination**.
11. Change the format to **FASTA**
12. Click on **Create File** to download the FASTA File.
13. **Repeat** this process for a gene interesting to you.

Sequence Retrieval

 **National Library of Medicine**
National Center for Biotechnology Information

Log in

Nucleotide Nucleotide Search

Advanced Help

GenBank Send to Annotations shown

Fusarium oxysporum strain 4287 Lcc9 (lcc9) gene, complete cds

GenBank: EF990898.1
[FASTA](#) [Graphics](#)



Go to: ▼

LOCUS	EF990898	3669 bp	DNA	linear	PLN 22-JUL-2007
DEFINITION	Fusarium oxysporum strain 4287 Lcc9 (lcc9) gene, complete cds.				
ACCESSION	EF990898				
VERSION	EF990898.1				
KEYWORDS	.				
SOURCE	Fusarium oxysporum				
ORGANISM	Fusarium oxysporum Eukaryota; Fungi; Dikarya; Ascomycota; Pezizomycotina; Sordariomycetes; Hypocreomycetidae; Hypocreales; Nectriaceae; Fusarium; Fusarium oxysporum species complex.				
REFERENCE	1 (bases 1 to 3669)				
AUTHORS	Cordoba C.D. and Roncero,M.I.G.				
TITLE	Isolation and characterization of laccases genes from the plant pathogenic fungus Fusarium oxysporum				
JOURNAL	Unpublished				
REFERENCE	2 (bases 1 to 3669)				
AUTHORS	Cordoba C.D. and Roncero,M.I.G.				
TITLE	Direct Submission				
JOURNAL	Submitted (26-JUN-2007) Genetic, University of Cordoba, Campus de Rabanales, Cordoba E-14071, Spain				
FEATURES	Location/Qualifiers				
source	1..3669				

☒ Complete Record
☐ Coding Sequences
☐ Gene Features

Choose Destination
☒ File ☐ Clipboard
☐ Collections ☐ Analysis Tool

Download 1 item.
Format
FASTA
☐ Show GI
Create File

Recent activity
Turn Off Clear
 Fusarium oxysporum strain 4287 Lcc9 (lcc9) gene, complete cds
 lcc9 (27)
Nucleotide


Abstract geometric shapes in the top-left corner, including a light blue triangle, a teal triangle, a green triangle, a red triangle, and a purple triangle, all pointing towards the center.

BLAST

Abstract geometric shapes in the bottom-right corner, consisting of several overlapping light gray triangles pointing towards the center.

BLAST



- BLAST stands for Basic Local Alignment Search Tool.
 - It identifies similarities between biological sequences by comparing nucleotide or protein sequences to a database of sequences.
- 



Important update

The **ClusteredNR** database is now the **default** Protein BLAST database. [Learn more about ClusteredNR](#)

Basic Local Alignment Search Tool

BLAST finds regions of similarity between biological sequences. The program compares nucleotide or protein sequences to sequence databases and calculates the statistical significance. [Learn more](#)

NEW



Mon, 21 Jul 2025

Here are a few highlights in our latest BLAST+ release:

Download BLAST+ 2.17.0 now!

[More BLAST news...](#)

Web BLAST



BLAST Genomes

Enter organism common name, scientific name, or tax id

Search

Activate
Go to Settings

BLAST



- **BLASTn (Nucleotide BLAST):** compares one or more nucleotide query sequences to a subject nucleotide sequence or a database of nucleotide sequences.
- **BLASTx (translated nucleotide sequence searched against protein sequences):** compares a nucleotide query sequence that is translated in six reading frames against a database of protein sequences.
- And **many others.**

BLAST



1. Choose **blastn**.
2. Paste the accession number **EF990898.1** in the Enter Query Sequence box.
3. Select **Nucleotide collection nr/nt** under Database in Choose Search Set section.
4. Under Program Selection choose **Highly similar sequences(megablast)**.
5. Set the max target sequences under General Parameters to **50**.
6. Click the **BLAST** button.



Important update

The **ClusteredNR** database is now the **default** Protein BLAST database. [Learn more about ClusteredNR](#)

[← Edit Search](#)[Save Search](#)[Search Summary ▾](#)[? How to read this report?](#)[▶ BLAST Help Videos](#)[↶ Back to Traditional Results Page](#)

i Your search is limited to records **include: entrez query: complete**

Job Title	EF990898:Fusarium oxysporum strain 4287 Lcc9...
RID	AW74CD1E014 Search expires on 08-28 04:43 am Download All ▾
Program	BLASTN ? Citation ▾
Database	nt See details ▾
Query ID	EF990898.1
Description	Fusarium oxysporum strain 4287 Lcc9 (lcc9) gene, comple ...
Molecule type	nucleic acid
Query Length	3669
Other reports	Distance tree of results MSA viewer ?

[Descriptions](#)[Graphic Summary](#)[Alignments](#)[Taxonomy](#)

Filter Results

Organism *only top 20 will appear* ☐ exclude

[+ Add organism](#)

Percent Identity

 to

E value

 to

Query Coverage

 to [Filter](#)[Reset](#)

Summary

Filter

Activate Windows
Go to Settings to activate Windows.

Feedback



Descriptions

Graphic Summary

Alignments

Taxonomy

Sequences producing significant alignments

Download

Select columns

Show

50

select all

5 sequences selected

GenBank

Graphics

Distance tree of results

MSA Viewer

	Description	Scientific Name	Max Score	Total Score	Query Cover	E value	Per. Ident	Acc. Len	Accession
<input checked="" type="checkbox"/>	Fusarium oxysporum strain 4287 Lcc9 (lcc9) gene, complete cds	Fusarium oxysporum	6761	6761	100%	0.0	100.00%	3669	EF990898.1
<input checked="" type="checkbox"/>	Fusarium acuminatum isolate 1A chromosome 5	Fusarium acuminatum	1760	2238	89%	0.0	79.57%	4969144	CP151264.1
<input checked="" type="checkbox"/>	Fusarium graminearum chromosome 1, complete genome	Fusarium graminearum	1249	1249	51%	0.0	78.90%	11760891	HG970332.2
<input checked="" type="checkbox"/>	Fusarium asiaticum strain KCTC 16664 chromosome 1	Fusarium asiaticum	1219	1219	51%	0.0	78.53%	12002067	CP088257.1
<input checked="" type="checkbox"/>	Colletotrichum jinshuiense isolate Esh8 chromosome 2	Colletotrichum jinshuiense	115	115	3%	3e-20	87.88%	6134399	CP150805.1

Activate Windows

Go to Settings to activate Windows.

Activate Windows
Go to Settings to activate



Descriptions

Graphic Summary

Alignments

Taxonomy


 hover to see the title  click to show alignments


Alignment Scores

 < 40


 40 - 50

 50 - 80


 80 - 200

 >= 200



5 sequences selected 

Distribution of the top 6 Blast Hits on 5 subject sequences





Download ▾ GenBank Graphics Sort by: E value ▾

▾ Next ▲ Previous ◀ Descriptions

Fusarium acuminatum isolate 1A chromosome 5

Sequence ID: [CP151264.1](#) Length: 4969144 Number of Matches: 2

Range 1: 2399923 to 2402423 [GenBank](#) [Graphics](#)

▾ [Next Match](#) ▲ [Previous Match](#)

Score	Expect	Identities	Gaps	Strand
1760 bits(953)	0.0	2013/2530(80%)	51/2530(2%)	Plus/Minus
Query 986	TTTGCTGCCCCGCGCGTTAGCTGCGACGGTGTCTTATGACTTTTCTATTGATTGGGTTTCG			1045
Sbjct 2402423	TTTGATGCCGCGCGCTTTGGCTGCTACGGTGTCTTATGATTTCACTATTGAATGGGTACG			2402364
Query 1046	AGCAAATCCAGATGGCGCGTTTGAGAGGTCGACGATAGGCATTAATAGAGAGTGGCCGAT			1105
Sbjct 2402363	AGCGAATCCTGATGGCGCCTTTGAGAGGCCTACGATTGGCATCAATGGGCGGTGGCCGAT			2402304
Query 1106	ACCGAGGATTGAAGCGAGTATTGGGGATACGGTTTTGGTTTATGTGAGGAATAATTTGGG			1165
Sbjct 2402303	TCCCAGGATCGAGGCGACTGTGGGTGATACGATTTTGGTGAATGCGAGGAATAATCTGGG			2402244
Query 1166	GAATCAGTCTACGAGTTTGCATTTTCATGGGCTTTTCATGAATGGCTCGAATCATATGGA			1225
Sbjct 2402243	GAATCAGTCCACGTCGTTGCATTTTCACGGTCTGTTTATGAATGGTTCAAACCATATGGA			2402184
Query 1226	TGGGCCGTCGCAGGTTACGCAGTGCCTATTCAACCTGGAGAGTCATTTCTCTATAACTT			1285
Sbjct 2402183	TGGGCCGTCGCAGGTTACGCAGTGCCTATTCAACCTGGAGAGTCATTTCTCTATAACTT			2402124

The image features a white background with decorative geometric elements. In the top-left corner, there is a cluster of overlapping triangles in shades of blue, teal, green, and purple. In the bottom-right corner, there is a cluster of overlapping triangles in various shades of gray.

Quality Control

Quality Control in Galaxy



Obtain **FASTQ** files:

1. Go to [Galaxy](#)
2. In the tool Bar, click on **Get Data**
3. Choose “**Faster Download and Extract Reads in FASTA/Q format from NCBI SRA**”.
4. In the Accession tab, write the accession number of the fastq file: **SRR18453616**.

Quality Control in Galaxy

The screenshot displays the Galaxy web interface. On the left is a sidebar with navigation icons for Upload, Tools, Workflows, Workflow Invocations, Interactive Tools, Visualization, Histories, History, Notifications, and More. The main panel is titled 'All Tools' and contains a search bar and a list of tools. The tool 'Faster Download and Extract Reads in FASTQ format from NCBI SRA' (Galaxy Version 3.1.1+galaxy1) is selected. The tool parameters section shows 'select input type' set to 'SRR accession' and 'Accession' set to 'SRR18453616'. A note states: 'Must start with SRR, DRR or ERR, e.g. SRR925743, ERR343809'. The 'Advanced Options' section is expanded, showing 'Email notification' set to 'No' and 'Attempt to re-use jobs with identical parameters?' set to 'No'. A 'Run Tool' button is visible at the bottom of the tool configuration panel.

Galaxy

Using 7% of 2

Faster Download and Extract Reads in FASTQ format from NCBI SRA
(Galaxy Version 3.1.1+galaxy1)

Tool Parameters

select input type

SRR accession

Accession *

SRR18453616

Must start with SRR, DRR or ERR, e.g. SRR925743, ERR343809

Advanced Options

Additional Options

Email notification

☐ No

Send an email notification when the job completes.

Attempt to re-use jobs with identical parameters?

☐ No

This may skip executing jobs that you have already run.

Run Tool

Help


Quality Control in Galaxy



Perform **QC**:

1. On the search bar, type **fastqc**.
2. **Choose** the desired fastq file (paired end) in the raw read tab.
3. **Leave** all other tabs unchanged.

Quality Control in Galaxy



The screenshot displays the Galaxy web interface for configuring the FastQC tool. The left sidebar contains navigation links: Upload, Tools, Workflows, Workflow Invocations, Interactive Tools, Visualization, Histories, History, and Notifications. The main panel shows the 'FastQC Read Quality reports (Galaxy Version 0.74+galaxy1)' tool configuration. The 'Tool Parameters' section includes:

- Raw read data from your current history ***: A dropdown menu showing selected datasets: '5: (hidden) SRR18453616:forward' and '6: SRR18453616:reverse'. Below the dropdown are links for 'accepted formats' and 'switch to column select'. A warning message states: 'This is a batch mode input field. Individual jobs will be triggered for each dataset.'
- Contaminant list - optional**: A dropdown menu currently set to 'Nothing selected'. Below it are links for 'accepted formats' and a description: 'tab delimited file with 2 columns: name and sequence. For example: Illumina Small RNA RT Primer CAAGCAGAAGACGGCATACGA'.
- Adapter list - optional**: A dropdown menu currently set to 'Nothing selected'. Below it are links for 'accepted formats' and a description: 'List of adapters adapter sequences which will be explicitly searched against the library. It should be a tab-delimited file with 2 columns: name and sequence. (--adapters)'.
- Submodule and Limit specifying file - optional**: A dropdown menu currently set to 'Nothing selected'. Below it are links for 'accepted formats' and a description: 'a file that specifies which submodules are to be executed (default=all) and also specifies the thresholds for the each submodules warning parameter'.
- Disable grouping of bases for reads >50bp**: A toggle switch currently set to 'No'.

The top right of the interface shows a status bar with 'Using 7% of 250'.

Quality Control in Galaxy



4. Once it runs, two files are generated, a raw data file and a Web Page file.
5. **View** the result by clicking on the web page file produced.

Quality Control Result

Summary

- ✓ [Basic Statistics](#)
- ✓ [Per base sequence quality](#)
- ✓ [Per sequence quality scores](#)
- ✗ [Per base sequence content](#)
- ✓ [Per sequence GC content](#)
- ✓ [Per base N content](#)
- ✓ [Sequence Length Distribution](#)
- ✗ [Sequence Duplication Levels](#)

✓ Per base sequence quality



Quality Control (MultiQC)



Why: It helps us to obtain a more intuitive comparison


1. On the search bar, type **multiqc**
2. On the “Which tool was used generate logs?” tab, choose **Fastqc**
3. Then click on “**Insert FastQC output**”

Quality Control (MultiQC)



4. Type of output is raw data
5. Add the raw data files generated earlier
6. **Leave** all other parameters at default
7. **Run** tool
8. **View** the result by clicking on the web page file produced

Quality Control (MultiQC)



The screenshot displays the MultiQC tool interface within the Galaxy web portal. The top navigation bar includes the Galaxy logo, a search bar, and a status indicator showing 'Using 13% of 256 MB'. The left sidebar contains navigation links for Upload, Tools, Workflows, Workflow Invocations, Interactive Tools, Visualization, Histories, History, and Notifications. The main content area is titled 'MultiQC' and includes a search bar with the text 'multiqc'. Below the search bar, there is a 'Show Sections' button. The 'Tool Parameters' section is currently empty. The 'Results' section is titled '1: Results' and contains a dropdown menu for 'Which tool was used generate logs?' with 'FastQC' selected. Below this, there is a 'FastQC output' section with a dropdown for 'Type of FastQC output?' set to 'Raw data'. The 'FastQC output' section also displays a list of outputs: '8: FastQC on data 5: RawData' and '10: FastQC on data 6: RawData'. At the bottom of the results section, there is a '+ Insert FastQC output' button and a '+ Insert Results' button. The footer of the interface shows 'Image content - optional'.

Quality Control Result

The screenshot displays the Galaxy web interface. On the left is a sidebar with navigation icons for Upload, Tools, Workflows, Workflow Invocations, Interactive Tools, Visualization, Histories, History, Notifications, and More. The main panel shows the '11: MultiQC on data 8 and data 10: Webpage' job, which is 4.7 MB in size and in HTML format. Below the job title are buttons for Preview, Visualize, Details, and Edit. The MultiQC tool version is v1.27. A left-hand menu lists various sections: General Stats, FastQC, Sequence Counts, Sequence Quality Histograms, Per Sequence Quality Scores, Per Base Sequence Content, Per Sequence GC Content, Per Base N Content, Sequence Length Distribution, Sequence Duplication Levels, Overrepresented sequences by sample, and Top overrepresented sequences. The main content area displays the MultiQC logo and a description: 'A modular tool to aggregate results from bioinformatics analyses across many samples into a single report.' It also shows the report generation date and time (2025-08-25, 03:54 UTC) and the data source path. A 'Welcome!' message with a 'Watch a tutorial video' button is visible. Below this is the 'General Statistics' section, which includes buttons for 'Copy table', 'Configure columns', 'Scatter plot', 'Violin plot', and 'Export as CSV...'. A table is partially visible with columns for 'Sample Name', 'Dups', 'GC', and 'Seqs'.

Galaxy Using 13% of 25

11: MultiQC on data 8 and data 10: Webpage ok

4.7 MB
format **html** database ? size **4.7 MB**

MultiQC aggregate results from bioinformatics analyses into a single report

multiqc v1.27

General Stats

FastQC

Sequence Counts

Sequence Quality Histograms

Per Sequence Quality Scores

Per Base Sequence Content

Per Sequence GC Content

Per Base N Content

Sequence Length Distribution

Sequence Duplication Levels

Overrepresented sequences by sample

Top overrepresented sequences

multiqc

A modular tool to aggregate results from bioinformatics analyses across many samples into a single report.

Report generated on 2025-08-25, 03:54 UTC based on data in: /corral4/main/jobs/070/149/70149916/working/multiqc_WDir

Welcome! Not sure where to start? don't show again

Watch a tutorial video (6:06)

General Statistics

Copy table Configure columns Scatter plot Violin plot Export as CSV...

Showing $2\frac{1}{2}$ rows and $3\frac{1}{6}$ columns.

Summarize table

Sample Name	Dups	GC	Seqs
-------------	------	----	------

https://usegalaxy.org/datasets/f9cad7b01a47213585b73fee9118e9d4

Why QC is Important

- Ensures data integrity before downstream analysis.
- Detects contamination, errors, or poor-quality sequences.
- Prevents misleading results.

Why QC is Important

- Ensures data integrity before downstream analysis.
- Detects contamination, errors, or poor-quality sequences.
- Prevents misleading results.